

Datasheet

Product Name: EEG Cap 32 channels based on CA-065 identical CW-0434
Product Code: CW-3484
Product Version: 1

Manufacturer:

 eemagine Medical Imaging Solutions
 Gubener Str. 47
 10243 Berlin
 Germany
 Phone +49 (0)30 2904 8404
 Fax +49 (0)30 2904 8405
 E-Mail support@eemagine.com
 Web www.eemagine.com

Classification: This custom work is not CE marked and strictly for research only.
Any classification that may apply to a reference product listed in this document for comparison and completeness of the custom work specification does NOT apply to this custom work!

Compatibility: Not specified. Pinning defined only.
Compatibility of custom work is neither verified nor validated. It is only implied by similarity to given reference products or by design.

Description: EEG cap, 32 electrodes, layout 10/20, sintered Ag/AgCl, 1x Redel connector, unshielded cable 1.8m, incl. 5 openings for HPI coils. All electrode labels comprise channel numbers and electrode names.

CW-3484 is designed to comply with the requirements to be used inside MEG environments.


Available sizes: C
(if applicable)

Specifications:**Cap Layout**

left

frontal

right

Openings for HPI coils 1-5

Pinning scheme:

Channel Number	Electrode	Pin Numbers
1	Fp1	1
2	Fp2	2
3	F7	3
4	F3	4
5	Fz	5
6	F4	6
7	F8	7
8	T7	8
9	C3	9
10	Cz	10
11	C4	11
12	T8	12
13	P7	13
14	P3	14
15	Pz	15
16	P4	16
17	P8	17
18	O1	18
19	O2	19
20	Fpz	20
21	FC5	21
22	FC1	22
23	FC2	23
24	FC6	24
25	CP5	25
26	CP1	26
27	CP2	27
28	CP6	28
29	M1	29
30	POz	30
31	M2	31
32	Oz	32

Cabling	
Cable length	1.8m
Cable	unshielded
Connector	
Connector type	Redel (32 pins)
Connector gender	Male

Warning:

Proper use of the CW-3484 depends on careful reading of all instructions including the descriptions and labels that come with or on the devices. Inaccurate measurements may be caused by incorrect use of the device. Non-compliance with warnings and safety regulations may result in severe personal injury and total loss of equipment.

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English: Information for waveguard™ EEG cap and accessories

Manufacturer: eemagine Medical Imaging Solutions GmbH, Gubener Str. 47, 10243 Berlin, Germany

The waveguard EEG cap is a medical device of CE class I according to EU regulation 93/42/EEC. In the USA, waveguard EEG cap is a class II device according to the CFR, FDA 510(k) No. K110223. In Canada, certain cap types are registered as medical device class II according to the MDR under MDL no. 88779.

For product registrations in other countries send your inquiry to sales@eemagine.com.

Safety information:

Intended use

- waveguard EEG caps are intended for recording of electroencephalographic (EEG) signals, including event-related potentials (ERP) and evoked potentials (EP) on subjects or patients with intact scalp in combination with medical-grade EEG recording equipment (CE class IIa certified inside the European Union).
- The cap must be dry and clean when used. The application of caps on subjects or patients should only be carried out by medically trained personnel (EEG technician). The technician must remain present during the entire recording session to prevent patient risks such as strangulation. Special infection control requirements may apply.

Contra-indications

- The caps should not be used on patients with skin injuries or problems such as blisters, burns, wounds from operation or other superficial wounds or transmissible diseases, such as – but not limited to – hepatitis, human immunodeficiency virus (HIV) or Creutzfeldt-Jacob disease.
- Do not use the cap in a Magnetic Resonance Imaging (MRI) scanner or in any other high-intensity electromagnetic field radiation environment unless the type of the cap expressly states the compatibility with such a device.
- In the event of defibrillation, disconnect the cap and/or electrodes from the EEG recording equipment first.

Side effects

- No side effects are known if properly applied. Please apply the correct size of the cap. Keep in mind that a too small cap may cause headache.
- Slight red marks on the exposed skin (e.g. at the forehead) due to contact pressure are a regular effect of wearing the cap and disappear within short time.

Important technical notes

- Sintered Ag/AgCl electrode elements do not need any rechlorification.
- Do not use incompatible electrode materials like tin, platinum or stainless steel in one recording.
- Do not use hard brushes or any other hard materials (in particular no metal tools such as screw drivers) to clean electrode pins.
- Apply always the correct size of cap (related to head circumference).
- Make sure that the cap is completely dry before use.
- Do not clean the cap using a supersonic cleaning device.
- Do not keep the cap submerged in water for prolonged periods of time (maximum one to two hours).
- Do not expose the cap to hot environments over 50° Celsius (or short rinsing at 60° Celsius). The cap fabrics as well as the wiring do not withstand higher temperatures. Hot-air drying may significantly reduce the lifetime of the cap.
- Do not use any autoclaving procedure on the cap.
- Make sure that the cap connectors always remain dry.
- Please dry the cap ideally in an upside down position, with the cable being placed higher than the wet cap, in order to prevent adapter and tubing to be contaminated with water.
- Always disconnect caps holding the connector housing. Never drag at the cable! Note that dragged-out cables are not covered by the cap warranty.

The User Guide contains additional information.

See www.ant-neuro.com/support, or contact your local distributor to obtain a copy.

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English: Safety information for waveguard™ Neonatal cap

Special warnings apply: this cap can be used for neonates inside an incubator. For a correct application, use the correct cap size! You must be extra careful with the skin preparation and gelling procedures, and may consider using extra-thick gel that allows pre-gelling the cap prior to placing the cap on the baby's head. The cap may require a special washing procedure.

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
Datasheet

Product Name: EEG cap, 64 ch, identical to CW-0435
Product Code: CW-3485
Product Version: 1


Manufacturer:
 eemagine Medical Imaging Solutions
Gubener Str. 47
10243 Berlin
Germany
Phone +49 (0)30 2904 8404
Fax +49 (0)30 2904 8405
E-Mail support@eemagine.com
Web www.eemagine.com

 **OQC passed**
for part(s) with SERIAL number as indicated on product / package

- Classification:** This custom work is not CE marked and strictly for research only.
Any classification that may apply to a reference product listed in this document for comparison and completeness of the custom work specification does NOT apply to this custom work!
- Compatibility:** Not specified. Pinning defined only.
Compatibility of custom work is neither verified nor validated. It is only implied by similarity to given reference products or by design.
- Description:** EEG cap, 64 electrodes, layout 10/10, sintered Ag/AgCl, 2x Redel connector, unshielded cable 1.8m, incl. 5 openings for HPI coils. All electrode labels comprise channel number and electrode name.
- CW-3485 is designed to comply with the requirements to be used inside MEG environments.
- Available sizes:** S, M
(if applicable)

Specifications:
Cap Layout

 Openings for HPI coils

Pinning scheme:

Connector 1 – yellow bend relief		
Channel Number	Electrode	Connector Pin
1	Fp1	1
2	Fp2	2
3	F7	3
4	F3	4
5	Fz	5
6	F4	6
7	F8	7
8	T7	8
9	C3	9
10	Cz	10
11	C4	11
12	T8	12
13	P7	13
14	P3	14
15	Pz	15
16	P4	16
17	P8	17
18	O1	18
19	O2	19
20	Fpz	20
21	FC5	21
22	FC1	22
23	FC2	23
24	FC6	24
25	CP5	25
26	CP1	26
27	CP2	27
28	CP6	28
29	FT9	29
30	POz	30
31	FT10	31
32	Oz	32

Connector 2 – green bend relief		
Channel Number	Electrode	Connector Pin
33	AF7	1
34	AF3	2
35	AFz	3
36	AF4	4
37	AF8	5
38	TP9	6
39	F1	7
40	F2	8
41	TP10	9
42	FC3	10
43	M1	11
44	FC4	12
45	C5	13
46	C1	14
47	C2	15
48	C6	16
49	CP3	17
50	M2	18
51	CP4	19
52	P5	20
53	P1	21
54	P2	22
55	P6	23
56	PO3	24
57	PO4	25
58	FT7	26
59	FT8	27
60	TP7	28
61	TP8	29
62	PO7	30
63	PO8	31
64	CPz	32

Cabling	
Cable length	1.8m
Cable	unshielded
Connector type	Redel (32 pins)
Connector gender	Male

Warning:

Proper use of the CW-3485 depends on careful reading of all instructions including the descriptions and labels that come with or on the devices. Inaccurate measurements may be caused by incorrect use of the device. Non-compliance with warnings and safety regulations may result in severe personal injury and total loss of equipment.

Disclaimer: We have attempted to write this document as accurately as possible. However, mistakes are bound to occur, and we reserve the right to make changes to the products, which may render parts of this document invalid. No part of this document may be copied or reproduced without the explicit permission of the authors.

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